

24TH
WFHSS
CONGRESS
BRUSSELS
18-21 OCTOBER
2023

Hygiene safety and resource savings in automated cleaning – A contradiction?

Name: Dennis Eisert
Ulf Klung
Matthias Tschoerner
Bastian Wulff

Affiliation: Chemsiche Fabrik Dr.
Weigert
GmbH & Co. KG

wfhss
World Federation for
Hospital Sterilisation Sciences

Sustainability

Hannß Carl von Carlowitz, (1645 - 1714)

Sylvicultura Oeconomica, 1713, pp. 105–106

"... like such a conservation and cultivation of wood / that there is a continuous, constant and sustainable use / because it is an indispensable thing / without which the country cannot remain in its being."

<https://digital.slub-dresden.de/werkansicht/df/85039/9>

Print: Leipzig: Braun, 1713, Online-Ed.: Dresden: SLUB, 2013

<https://en.wikipedia.org/wiki/Sustainability>



Sustainability



... means meeting the needs of the present in a way that does not limit the opportunities of future generations.

Federal Ministry for Economic Cooperation and Development

<https://www.bmz.de/de/service/lexikon/nachhaltigkeit-nachhaltige-entwicklung-14700>

Energy and CO₂ balances

Examples of energy efficient hospitals

| Energy efficiency standard | Consumption characteristics of electrical primary energy arithmetic . Middle [kWh/m ² a] | Electrical energy consumption figures GHG emissions arithm . Middle [kg/m ² a] |
|---------------------------------------|---|---|
| Administration building | 77 | 20 |
| Primary care hospitals (-250 BE) | 171 | 45 |
| Hospitals standard care (251-450 BE) | 217 | 57 |
| Hospitals maximum care (from 1000 BE) | 134 | 35 |

Energy consumption characteristics of energetically highly efficient buildings, final report, KWEFF2015, C. Zeine, S. Gausling , C. Cassebaum , M. Gebhardt, N. Goldau, J. Peters, Society for Energy Planning and System Analysis mbH , Münster, September 2015

Environmental aspects in the SPD/CSSD

Latent questions for more than 10 years (Selection)

Reducing the environmental footprint of gastrointestinal endoscopy: European Society of Gastrointestinal Endoscopy (ESGE) and European Society of Gastroenterology and Endoscopy Nurses and Associates (ESGENA) Position Statement; de Santiago Enrique Rodríguez et al. Reducing the... Endoscopy 2022; 54:797-826

wfhss 2015 Lille: Kristian Hemström, study on the impact from production to disposal on energy consumption and the climate: **The use of single-use instruments releases 600 times as many greenhouse gases (GHG) consumes 17 times as much energy as the use of reusable instruments;** Report: G. Westermann, Central Sterilization, 2015, 388

Eco-Sterilization: A challenge for all of us; Christophe Lambert, 7th Swiss conference SGSV, 2011

Sustainable development in processing, 32nd French National Sterilization Days in Lille, April 28th and 29th, 2010; Report: G. Westermann, Central Sterilization, 2010, 144

Material and energy balance

Model CSSD – 4 WD, 1 Trolley washer, 225 operating days, 1 shift operation

| | Chemistry ml/l | Electricity cycle kWh | water cycle L | Cycles Year | Chemistry Year L | Electricity Year kWh | Water Year L |
|----------------------|-------------------|-----------------------------|---------------------|----------------|------------------------|--------------------------------|--------------------|
| WD | 5 | 10 | 120 | 7200 | 4320 | 72000 | 864000 |
| Trolley washer | 3 | 13.3 | 160 | 2700 | 1296 | 36000 | 432000 |
| Total consumption | | | | | 5616 | 108000 | 1296000 |
| Costs | | | | | | 37800 € | 51840 € |
| GHG emissions | | | | | | 54540 kg CO₂ | |

Operating days: 225 p. year, WD: 8 cycles p. day, Trolley washer 12 cycles p. day


Electricity costs € 0.35/kWh, deionized water costs € 0.04/L

Greenhouse gases CO₂ equivalent: 0.505 [kg/kWh], IWU 2020, Institute for Housing and Environment, 26-02-2020

Energy saving through optimized program sequences



WD Standard program sequence



WD Standard program sequence


| | |
|--|---|
| WD: | MMM UniClean PL II |
| Test spicemen: | Crile clamp [1] |
| Test soil: | Coagulated blood [1, 2] |
| Pre-rinse: | 2 min, cold |
| Cleaning: | 10 min at 55 °C 6 ml/l mild-alkaline enzymatic cleaner |
| Neutralization: | 2 min, cold |
| Intermediate rinse, Thermal disinfection, Drying | |

| Parameter | Result per cycle |
|--|--------------------|
| Water | 178 L |
| Time | 49 mins |
| Electricity | 10 kWh |
| Process chemistry | 220 ml |
| Cleaning performance (residual protein) | Ø 14.4 µg ± 5.9 µg |

[1] Guideline compiled by DGKH, DGSV and AKI for the validation and routine monitoring of automated cleaning and thermal disinfection processes for medical devices; 5th Edition. Central Service 2017; 25, Suppl.
 [2] ISO 15883-5:2021 Table A.1, 1st Row; [3] AAMI ST98 (2023)



WD energy optimized program sequence



WD energy optimized program sequence

| | |
|----------------|--|
| WD: | MMM UniClean PL II, electrical heated |
| Test specimen: | Crile clamp [1] |
| Test soil: | Coagulated blood [1, 2] |
| Pre-rinse: | 2 min, cold, 1 ml/l optimized detergent |
| Cleaning: | 10 min at 35 °C |
| | 1 ml/l optimized detergent, dos. 25 °C |
| | Intermediate rinse, Thermal disinfection, Drying |


| Parameter | Result per cycle |
|---|------------------------------|
| Water | 141 L |
| Time | 36 mins |
| Electricity | 8 kWh |
| Process chemistry | 75 ml |
| Cleaning performance (residual protein) | Ø 15.6 µg ± 11.5 µg * |

* < 3,0 µg/cm²


[1] Guideline compiled by DGKH, DGSV and AKI for the validation and routine monitoring of automated cleaning and thermal disinfection processes for medical devices; 5th Edition. Central Service 2017; 25, Suppl.
[2] ISO 15883-5:2021 Table A.1, 1st Row; [3] AAMI ST98 (2023)



WD energy optimized program sequence



Cleaning performance Energy optimized



➔ **Consistent cleaning result despite energy savings**

Real use - WD Validation France

WD: Belimed WD 290, electrical heated

Load (selection):

- Surgical instruments

- MIS instruments

- Gynecological instruments

- Ophthalmic instruments

- Anesthetic equipment

- Rigid endoscopes / optics

- Container (anodised aluminium)

- DaVinci Optics



Real use - WD Validation France

WD: Belimed WD 290, electrical heated

Contamination from real use (selection):

Blood, Protein

Denatured blood (e.g. RF-Technology)

Mucus, blood-mucus mixtures (e.g. gynecology)

Mucus (ENT)

Dental dirt, mouth, jaw and palate

Bone meal, dentin

Contamination from ophthalmology _

Antiseptic Residues (containing iodine)




Real use - WD Validation France

| | Commonly used mild-alkaline enzymatic detergent | Optimized mild-alkaline enzymatic detergent |
|--------------------------|--|--|
| Pre-rinse | Unchanged | Unchanged |
| Cleaning | 4 ml/l, 10 min, 55 °C, dosing at 38 °C | 1 ml/l, 5 min, 45 °C, dosing at 25 °C |
| Intermediate rinse | Unchanged | Unchanged |
| Water | 35 l in the cleaning step | 35 l in the cleaning step |
| Process chemistry | 140ml | 35ml |


Real use - WD Validation France



Real use - WD Validation France



Cleaning performance Energy optimized



Consistent cleaning result despite energy savings

Real use - WD Validation – Service provider Germany

WD: Steelco DS 1000, steam heated

Load (selection):

- Surgical instruments
- MIS instruments
- Gynecological instruments
- Endoprosthetic sets
- Rigid endoscopes / optics
- Container
- Stationware (instruments, utensils)



Real use - WD Validation – Service provider Germany

WD: Steelco DS 1000, steam heated

Contamination from real use (selection):

Blood, Protein

Denatured blood (e.g. HF-technology)

Mucus, blood-mucus mixtures (e.g. gynecology)


Mucus (ENT)

Bone meal, dentin



Real use - WD Validation – Service provider Germany

| | Commonly used mild-alkaline enzymatic detergent | Optimized mild-alkaline enzymatic detergent |
|--------------------------|---|--|
| Pre-rinse | 2 min, cold | 2 min, cold, Dos. 0.5 ml/l |
| Cleaning | 6 ml/l, 10 min, 55 °C, dosing at 38 °C | 2 ml/l, 5 min, 35 °C, dosing at 20 °C |
| Intermediate rinse | Unchanged | Unchanged |
| Water | 170 l | 170 l |
| Process chemistry | 270 ml | 113 ml (23 ml + 90 ml) |



Real use - WD Validation – Service provider Germany

energy demand heated by steam,
calculated


$$\Delta Q = c * m * \Delta T$$

$$\Delta Q = 4,19 \frac{kJ}{kg * K} * 40 kg * 20 K = 3352 kJ$$

$$\Delta Q = 3352 kJ = 0,93 kWh$$


$$m = \frac{\Delta Q}{q_v} = \frac{3352 kJ}{2257 \frac{kJ}{kg}} = 1,48 kg \text{ Steam}$$

other energy savings per cycle,
measured




$$\Delta E = 0,033 kWh$$

Real use - WD Validation – Service provider Germany





Cleaning performance Energy optimized



Consistent cleaning result despite energy savings

Material and energy balance savings

Model CSSD – 4 WD, 1 Trolley washer, 225 operating days, 1 shift operation



| Saving | | | | |
|----------------|---------------------------------------|--------------------------------|------------------------------|--------------------------------|
| | Process chemicals p. cycle ml/l | Electricity p. cycle kWh | Process chemicals Year | Electricity Year |
| WD | 2.5 | 2 | 2160 L | 14400 kWh |
| Trolley washer | 2 | 2.7 | 864 L | 7200 kWh |
| Total savings | | | 3024 L | 21600 kWh |
| Costs | | | | € 7560 |
| GHG emissions | | | | 10908 kg CO₂ |

Operating days: 225 pa, WD: 8 cycles p. day – 3200 pa., Trolley washer: 12 batches pd . – 2700 pa
Electricity costs € 0.35/kWh, deionized water costs € 0.04/L

Greenhouse gases CO₂ equivalent: 0.505 [kg/kWh], IWU 2020, Institute for Housing and Environment, 26-02-2020

Other beneficial effects

- ✓ Less packaging material
- ✓ Lower transport costs and transport emissions
- ✓ Less storage space for replacement containers
- ✓ Fewer container replacement
- ✓ Less waste

Questions?



Thanks!